Nouryon

Statement

Efficacy of ARQUAD[®] MCB-50 versus influenza virus H1N1 (Swine Flu, Mexican flu)

The efficacy of ARQUAD[®] MCB-50 was tested in a suspension test following to the European standard DIN EN 14476:2005-08 under clean conditions (0,03% BSA).

ARQUAD[®] MCB-50 was tested as a 750 ppm, 500 ppm, and 250 ppm solution. The exposure times were 1 min, 5 min, 10 min and 15 min.

Conclusion: 750 ppm concentration of ARQUAD[®] MCB-50 (is 375 ppm active substance) is effective against influenza virus A/H1N1/X-179A at room temperature under clean conditions within an application time of 10 minutes.

SE-444 85 Stenungsund Sweden T +46 303 85000 www.nouryon.com